Liquichek™ Blood Gas Control Levels 1, 2 and 3

REF	571	Level 1	30 x 1.7 mL
	572	Level 2	30 x 1.7 mL
	573	Level 3	30 x 1.7 mL
	570X	Trilevel MiniPak	9 x 1.7 mL

EXP 2027-09-30

28780

Level 1	28781
Level 2	28782
Level 3	28783

ENGLISH

INTENDED USE

Liquichek Blood Gas Control is an assayed quality control intended for use in monitoring the precision of an individual laboratory's measurement of pH, pCO₂ and pO₂ by blood gas instrumentation.

SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Three levels of control are available to allow performance monitoring within the clinical range.

For customers in Germany: Quality control materials are required for assessment of laboratory performance as described in the "Guideline for Quality Assurance of Medical Laboratory Examinations following the German Medical Association" (Rili-BÄK regulation).

REAGENT

This product is a buffered bicarbonate solution in equilibrium with predetermined levels of oxygen, carbon dioxide and nitrogen.

Safety Data Sheet (SDS) available for professional users on www.bio-rad.com.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at 18 to 25°C. Once the control is opened, all analytes should be assayed immediately; discard the remaining material.

Avoid storage above 30°C. Freezing this control will not adversely affect performance.

This product is shipped under ambient conditions.

PROCEDURE

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit, or reagent being used.

For optimum performance, allow a minimum of 4 hours for the ampules to equilibrate to the temperature of 20 to 23°C. Before sampling, hold the ampule between the thumb and forefinger and shake the control for 10 seconds to completely mix the contents. Restore the control to the bottom of the ampule by gentle tapping. Protect fingers with a tissue or gloves to avoid injury. Snap the ampule open and immediately introduce the control into the instrument using one of the following techniques: direct aspiration, syringe transfer or capillary mode transfer.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

LIMITATIONS

1. This product should not be used past the expiration date.
2. This product is not intended for use as a standard.
3. This product is sensitive to most instrument related factors which cause unexpected analytical deviations. However, this control does not contain red blood cells and may not detect conditions that cause lysis of red blood cells, resulting in erroneous recovery of blood gas values.

ASSIGNMENT OF VALUES

The means and acceptable ranges reported in this insert were derived from replicate analyses and are specific for this lot of product. These tests were performed on several blood gas analyzers. For some instruments the assigned values are based on comparative studies for each make and model. Individual laboratory means should fall within the corresponding acceptable range, however, other factors such as procedural or reagent modifications may alter the means. Each laboratory should establish its own means and acceptable ranges and use the values provided only as guides.

SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a stabilized liquid product manufactured under rigid quality control standards. To obtain consistent ampule-to-ampule assay values, the control requires proper storage and handling as described.

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ENGLISH

▲ Data not available at the time of printing. Please inquire.

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Instrumentation Laboratory (IL), Lexington, Massachusetts
L. Eschweiler and Company, Kiel, Germany
Nova Biomedical, Waltham, Massachusetts

Radiometer of America, Westlake, Ohio
Roche Diagnostics Corporation, Indianapolis, Indiana
Siemens Healthcare Diagnostics Inc., Tarrytown, New York

Bio-Rad Laboratories comprehensive line of quality controls and QC data management solutions.

Autoimmune Controls
Blood Gas Controls
Cardiac Assessment Controls
Chemistry Controls
Coagulation Controls
Congenital/Pediatric Disease Controls

Diabetes/Hemoglobin Controls
Hematology Controls
Hepatitis & Retrovirus Controls
Immunoassay Controls
Immunology/Protein Controls
Molecular Controls
Sexually Transmitted Disease Controls

Specialty Infectious Disease Controls
Therapeutic Drug Monitoring Controls
Toxicology: Drugs-of-Abuse Controls
Toxicology: Specialty Controls
Urinalysis Controls
QC Data Management Solutions
External Quality Assurance Services (EQAS)

BIO-RAD

USA: 1-800-224-6723 | www.bio-rad.com/qualitycontrol

04/14

Mfd. For

UNITED STATES, Bio-Rad Laboratories
9500 Jeronimo Road, Irvine, CA 92618

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