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### 1.1、Experiment Properties 设置

按照下图所示选择：Genotyping → Taqman® Reagents → Standard → Include: Pre-PCR、Read Amplification、Post-PCR Read

What type of experiment do you want to set up?

Quantitation - Standard Curve    Quantitation - Relative Standard Curve    Quantitation - Comparative Ct ( $\Delta\Delta C_t$ )

Melt Curve    **✓ Genotyping**    Presence/Absence

Detect single nucleotide polymorphism variants of a target nucleic acid sequence in samples.

---

Which reagents do you want to use to detect the target sequence?

**✓ TaqMan® Reagents**    Other

The PCR reactions contain primers designed to amplify the target sequence and a TaqMan® probe designed to detect amplification of the target sequence.

---

Which ramp speed do you want to use in the instrument run?

**✓ Standard (~ 2 hours to complete a run)**

For optimal results with the standard ramp speed, Applied Biosystems recommends using standard reagents for your PCR reactions.

---

What do you want to include in the instrument run?

Include:  Pre-PCR Read     Amplification     Post-PCR Read

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Assign SNP Assay(s) to the Selected Wells.

Create New SNP Assay    Add Saved SNP Assay    Edit ▼

Assign	SNP Assay
<input type="checkbox"/>	SNP Assay 1

Save SNP Assay  
Edit SNP Assay...  
Save Selected SNP Assay  
Delete SNP Assay

Edit SNP Assay

Make changes below, then click "OK" to save your changes to the library. Click "Reset Fields" to undo all your changes.    **\*= Required**

SNP Assay Name: rs671    Color: ■    Assay ID:

Allele 1 Name or Base(s): G    Color: ■    Reporter: VIC    Quencher: NFQ-MGB

Allele 2 Name or Base(s): A    Color: ■    Reporter: FAM    Quencher: NFQ-MGB

Comments:

Reset Fields    OK    Cancel

按照下图所示设置 Assign Sample to the Selected Wells 中样本信息

**Assign Sample to the Selected Wells.**

Add New Sample   Add Saved Sample   Save Sample   Delete Sample

Assign	Sample	Color
<input type="checkbox"/>	gDNA 1	
<input type="checkbox"/>	Sample 2	
<input checked="" type="checkbox"/>	Blood 3	

按照下图所示设置 Select the dye to use as the passive reference 中 ROX 信息

**Select the dye to use as the passive reference.**

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**Assign SNP Assay(s) to the Selected Wells.**

Create New SNP Assay   Add Saved SNP Assay   Edit

Assign	SNP Ass...	Allele 1/Allele 2 Reporter	Task
<input checked="" type="checkbox"/>	rs671	VIC/FAM	Unknown

**Assign Sample to the Selected Wells.**

Add New Sample   Add Saved Sample   Save Sample   Delete Sample

Assign	Sample	Color
<input type="checkbox"/>	gDNA 1	
<input type="checkbox"/>	Sample 2	
<input checked="" type="checkbox"/>	Blood 3	

**Select the dye to use as the passive reference.**

ROX

**View Plate Layout**   View Well Table


Select Wells With: - Select Item -   - Select Item -

Show in Wells   View Legend

	1	2	3	4	5	6	7	8	9
A									
B				gDNA 1 U rs...					
C				Sample 2 U rs...					
D				Blood 3 U rs...					
E									
F									
G									
H									

Wells: U 3 Unknown   N 0 Negative Control   P 0 Positive Control

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
Reaction Volume Per  Well

Holding Stage: 95.0°C 10:00min

Cycling Stage: 95.0°C 00:10min 60.0°C 00:30min (Collect Data) Number of Cycles 40

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
Instrument Type: 7500 Fast Real-Time PCR System

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