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1.1、Experiment Properties 设置

按照下图所示选择：Genotyping → Taqman® Reagents → Standard → Include: Pre-PCR、Read Amplification、Post-PCR Read

What type of experiment do you want to set up?

Detect single nucleotide polymorphism variants of a target nucleic acid sequence in samples.

Which reagents do you want to use to detect the target sequence?

The PCR reactions contain primers designed to amplify the target sequence and a TaqMan® probe designed to detect amplification of the target sequence.

Which ramp speed do you want to use in the instrument run?

For optimal results with the standard ramp speed, Applied Biosystems recommends using standard reagents for your PCR reactions.

What do you want to include in the instrument run?

Include: Pre-PCR Read Amplification Post-PCR Read

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Assign SNP Assay(s) to the Selected Wells.

| Assign | SNP Assay |
|--------------------------|-------------|
| <input type="checkbox"/> | SNP Assay 1 |

Edit SNP Assay

Make changes below, then click "OK" to save your changes to the library. Click "Reset Fields" to undo all your changes. *= Required

SNP Assay Name:
 Color:
 Assay ID:

Allele 1 Name or Base(s):
 Color:
 Reporter:
 Quencher:

Allele 2 Name or Base(s):
 Color:
 Reporter:
 Quencher:

Comments:

按照下图所示设置 Assign Sample to the Selected Wells 中样本信息

Assign Sample to the Selected Wells.

Add New Sample Add Saved Sample Save Sample Delete Sample

| Assign | Sample | Color |
|-------------------------------------|----------|-------|
| <input type="checkbox"/> | gDNA 1 | |
| <input type="checkbox"/> | Sample 2 | |
| <input checked="" type="checkbox"/> | Blood 3 | |

按照下图所示设置 Select the dye to use as the passive reference 中 ROX 信息

Select the dye to use as the passive reference.

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Assign SNP Assay(s) to the Selected Wells.

Create New SNP Assay Add Saved SNP Assay Edit

| Assign | SNP Ass... | Allele 1/Allele 2 Reporter | Task |
|-------------------------------------|------------|----------------------------|---------|
| <input checked="" type="checkbox"/> | rs671 | VIC/FAM | Unknown |

Assign Sample to the Selected Wells.

Add New Sample Add Saved Sample Save Sample Delete Sample

| Assign | Sample | Color |
|-------------------------------------|----------|-------|
| <input type="checkbox"/> | gDNA 1 | |
| <input type="checkbox"/> | Sample 2 | |
| <input checked="" type="checkbox"/> | Blood 3 | |

Select the dye to use as the passive reference.

ROX

View Plate Layout View Well Table


Select Wells With: - Select Item - - Select Item -

Show in Wells View Legend

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 |
|---|---|---|---|---------------------|---|---|---|---|---|
| A | | | | | | | | | |
| B | | | | gDNA 1 U rs... | | | | | |
| C | | | | Sample 2 U rs... | | | | | |
| D | | | | Blood 3 U rs... | | | | | |
| E | | | | | | | | | |
| F | | | | | | | | | |
| G | | | | | | | | | |
| H | | | | | | | | | |

Wells: U 3 Unknown N 0 Negative Control P 0 Positive Control

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
Reaction Volume Per Well

Holding Stage: 95.0°C 10:00min

Cycling Stage: 95.0°C 00:10min 60.0°C 00:30min (Collect Data) Number of Cycles 40

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

Instrument Type: 7500 Fast Real-Time PCR System

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