

【英文名称】: COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)

【中文名称】: 凝血质控品（水平 2）

【货号】 (CAT. NUMBER): CG5022

【批号】 (LOT NUMBER): 405CG

【包装规格】 (SIZE): 12 x 1ml

【到期】 (EXPIRY): 2025-09-28

本产品的准备和使用详见试剂盒内说明书。

Full instructions for the use and preparation of this product are included with the kit.

质控品靶值和范围如下所示。

Target values and ranges are listed below.

ZH: 07 May '24 me

Method 所有仪器均值		COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)					凝血质控品 (水平2)		
批号 (Lot. No): 405CG		货号 (Cat. No): CG5022		到期 (Expiry): 2025-09-28					
包装规格 (Size): 12 x 1 ml		Range			范围				
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	1SD	2SD	Method 方法	
Activated Partial Thromboplastin Time (APTT)	活化部分凝血酶时间 (APTT)	seconds	47.9	38.3	57.5	4.80	9.60	Beijing Succeder aPTT	
		seconds	46.4	37.1	55.7	4.65	9.30	Bio-Ksel System APTTs	
		seconds	53.3	42.6	64.0	5.35	10.7	Biolabo Bio-CK APTT	
		seconds	73.5	58.8	88.2	7.35	14.7	Cobas aPTT MedS	
		seconds	57.9	46.3	69.5	5.80	11.6	DG Lab APTT	
		seconds	49.3	39.4	59.2	4.95	9.90	Diagon Dia-PTT Liquid	
		seconds	59.2	47.4	71.0	5.90	11.8	Diagon Dia-PTT Lyophilised	
		seconds	47.9	38.3	57.5	4.80	9.60	DIRUI aPTT	
		seconds	57.2	45.8	68.6	5.70	11.4	Erba Actime APTT	
		seconds	45.7	36.6	54.8	4.55	9.10	Fortress Diagnostics aPTT	
		seconds	57.8	46.2	69.4	5.80	11.6	Helena aPTT Si L Minus	
		seconds	64.8	51.8	77.8	6.50	13.0	HemosIL aPTT Lyophilised Silica	
		seconds	41.4	33.1	49.7	4.15	8.30	HemosIL APTT-SP liquid	
		seconds	61.1	48.9	73.3	6.10	12.2	HemosIL Synthasil	
		seconds	48.8	39.0	58.6	4.90	9.80	Human Hemostat APTT	
		seconds	54.0	43.2	64.8	5.40	10.8	Labitec aPTT	
		seconds	50.1	40.1	60.1	5.00	10.0	Media IVD aPTT	
		seconds	48.5	38.8	58.2	4.85	9.70	Mindray/Longisland APTT	
		seconds	48.4	38.7	58.1	4.85	9.70	MTI Diagnostics APTT	
		seconds	38.0	30.4	45.6	3.80	7.60	Renam APTT	
		seconds	62.0	49.6	74.4	6.20	12.4	Roche Cobas APTT	
		seconds	64.6	51.7	77.5	6.45	12.9	Roche Cobas APTT Screen	
		seconds	58.9	47.1	70.7	5.90	11.8	Sclavo APTT-S	
		seconds	46.2	37.0	55.4	4.60	9.20	Siemens/Dade Actin	
		seconds	47.5	38.0	57.0	4.75	9.50	Siemens/Dade Actin Cephaloplastin	
		seconds	46.2	37.0	55.4	4.60	9.20	Siemens/Dade Actin FS	
		seconds	41.2	33.0	49.4	4.10	8.20	Siemens/Dade Actin FSL	
		seconds	68.2	54.6	81.8	6.80	13.6	Siemens/Dade Actin SL	
		seconds	74.4	59.5	89.3	7.45	14.9	Siemens/Dade Pathromtin SL	
		seconds	59.5	47.6	71.4	5.95	11.9	Spinreact APTT	
		seconds	57.5	46.0	69.0	5.75	11.5	Stago aPTT	
		seconds	55.0	44.0	66.0	5.50	11.0	Stago Cephascreen	
		seconds	54.8	43.8	65.8	5.50	11.0	Stago CK Prest	
		seconds	56.8	45.4	68.2	5.70	11.4	Stago PTT Automate	
		seconds	58.6	46.9	70.3	5.85	11.7	Tcoag TriniCLOT aPTT S	
		seconds	46.6	37.3	55.9	4.65	9.30	Technoclone Siron aPTT	
		seconds	61.4	49.1	73.7	6.15	12.3	Teco TEClot APTT	
		seconds	50.8	40.6	61.0	5.10	10.2	Thalys APTT-EA	

Method 所有仪器均值		COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)					凝血质控品 (水平2)		
批号 (Lot. No): 405CG		货号 (Cat. No): CG5022		到期 (Expiry): 2025-09-28					
包装规格 (Size): 12 x 1 ml		Range			范围				
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	ISD	2SD	Method 方法	
Activated Partial Thromboplastin Time (APTT)	活化部分凝血酶时间 (APTT)	seconds	40.3	32.2	48.4	4.05	8.10	Wiener Lab APTT	
		seconds	45.5	36.4	54.6	4.55	9.10	Yumizen G APTT Liq	
Activated Partial Thromboplastin Time (APTT) Ratio	活化部分凝血酶时间 (APTT) 比值	RATIO	1.34	1.07	1.61	0.135	0.270	Bio-Ksel System APTTs	
		RATIO	1.81	1.45	2.17	0.180	0.360	Biolabo Bio-CK APTT	
		RATIO	1.72	1.38	2.06	0.170	0.340	Diagon Dia-PTT Lyophilised	
		RATIO	1.92	1.54	2.30	0.190	0.380	Erba Actime APTT	
		RATIO	2.09	1.67	2.51	0.210	0.420	Helena aPTT Si L Minus	
		RATIO	1.38	1.10	1.66	0.140	0.280	HemosIL APTT-SP liquid	
		RATIO	2.00	1.60	2.40	0.200	0.400	HemosIL Synthasil	
		RATIO	1.56	1.25	1.87	0.155	0.310	Human Hemostat APTT	
		RATIO	1.91	1.53	2.29	0.190	0.380	Labitec aPTT	
		RATIO	1.66	1.33	1.99	0.165	0.330	Media IVD aPTT	
		RATIO	1.63	1.30	1.96	0.165	0.330	MTI Diagnostics APTT	
		RATIO	2.14	1.71	2.57	0.215	0.430	Sciago APTT-S	
		RATIO	1.85	1.48	2.22	0.185	0.370	Siemens/Dade Actin FS	
		RATIO	1.54	1.23	1.85	0.155	0.310	Siemens/Dade Actin FSL	
		RATIO	2.28	1.82	2.74	0.230	0.460	Siemens/Dade Pathromtin SL	
		RATIO	1.79	1.43	2.15	0.180	0.360	Stago aPTT	
		RATIO	1.78	1.42	2.14	0.180	0.360	Stago Cephascreen	
		RATIO	1.78	1.42	2.14	0.180	0.360	Stago CK Prest	
		RATIO	1.66	1.33	1.99	0.165	0.330	Stago PTT Automate	
		RATIO	2.00	1.60	2.40	0.200	0.400	Tcoag TriniCLOT aPTT S	
		RATIO	2.31	1.85	2.77	0.230	0.460	Teco TEClot APTT	
		RATIO	1.90	1.52	2.28	0.190	0.380	Thalys APTT-EA	
		RATIO	1.47	1.18	1.76	0.145	0.290	Yumizen G APTT Liq	
Antithrombin III activity	抗凝血酶III活性	%	66.6	53.3	79.9	6.65	13.3	HemosIL Liquid Antithrombin	
		%	72.7	58.2	87.2	7.25	14.5	Roche Antithrombin	
		%	70.1	56.1	84.1	7.00	14.0	Roche Chromogenic Thrombin	
		%	68.5	54.8	82.2	6.85	13.7	Sciago ATIII	
		%	72.5	58.0	87.0	7.25	14.5	Siemens Innovance Antithrombin III	
		%	73.6	58.9	88.3	7.35	14.7	Siemens/Dade Antithrombin III	
		%	71.1	56.9	85.3	7.10	14.2	Siemens/Dade Berichrom ATIII	
		%	75.4	60.3	90.5	7.55	15.1	STACHrom ATIII	
		%	73.3	55.0	91.6	9.15	18.3	Thalys Antithrombin III	
		%	66.0	52.8	79.2	6.60	13.2	HemosIL Factor II deficient plasma	
Factor II activity	II因子活性	%	64.9	51.9	77.9	6.50	13.0	HemosIL RecombiPlasTin 2G	
		%	77.9	62.3	93.5	7.80	15.6	Siemens Factor II Deficient Plasma	

Method 所有仪器均值		COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)					凝血质控品 (水平2)		
批号 (Lot. No): 405CG		货号 (Cat. No): CG5022		到期 (Expiry): 2025-09-28					
包装规格 (Size): 12 x 1 ml		Range			范围				
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	1SD	2SD	Method 方法	
Factor II activity	II因子活性	%	74.7	59.8	89.6	7.45	14.9	Siemens/Dade Innovin	
		%	83.5	66.8	100	8.25	16.5	Siemens/Dade Thromborel S	
		%	65.6	52.5	78.7	6.55	13.1	Stago STA-NeoPTimal	
Factor IX activity	IX因子活性	%	83.1	66.5	99.7	8.30	16.6	HemosIL APTT-SP liquid	
		%	61.3	49.0	73.6	6.15	12.3	HemosIL Factor IX deficient plasma	
		%	65.6	52.5	78.7	6.55	13.1	HemosIL Synthasil	
		%	64.1	51.3	76.9	6.40	12.8	Siemens Factor IX Deficient Plasma	
		%	63.5	50.8	76.2	6.35	12.7	Siemens/Dade Actin FS	
		%	65.6	52.5	78.7	6.55	13.1	Siemens/Dade Actin FSL	
		%	70.2	56.2	84.2	7.00	14.0	Siemens/Dade Pathromtin SL	
		%	66.1	52.9	79.3	6.60	13.2	Stago Deficient / ImmunoDef IX	
Factor V activity	V因子活性	%	52.7	42.2	63.2	5.25	10.5	HemosIL Factor V deficient plasma	
		%	55.7	44.6	66.8	5.55	11.1	HemosIL RecombiPlasTin 2G	
		%	61.5	49.2	73.8	6.15	12.3	Siemens Factor V Deficient Plasma	
		%	50.1	40.1	60.1	5.00	10.0	Siemens/Dade Innovin	
		%	70.7	56.6	84.8	7.05	14.1	Siemens/Dade Thromborel S	
Factor VII activity	VII因子活性	%	54.2	43.4	65.0	5.40	10.8	Stago STA-NeoPTimal	
		%	94.4	75.5	113	9.30	18.6	HemosIL Factor VII deficient plasma	
		%	93.1	74.5	112	9.45	18.9	HemosIL RecombiPlasTin 2G	
		%	92.1	73.7	111	9.45	18.9	Siemens Factor VII Deficient Plasma	
		%	89.4	71.5	107	8.80	17.6	Siemens/Dade Innovin	
		%	95.2	76.2	114	9.40	18.8	Siemens/Dade Thromborel S	
Factor VIII activity	VIII因子活性	%	82.7	66.2	99.2	8.25	16.5	Stago Deficient Factor VII	
		%	41.3	33.0	49.6	4.15	8.30	HemosIL Factor VIII deficient plasma	
		%	43.7	35.0	52.4	4.35	8.70	HemosIL Synthasil	
		%	42.4	33.9	50.9	4.25	8.50	Siemens Factor VIII Deficient Plasma	
		%	41.8	33.4	50.2	4.20	8.40	Siemens/Dade Actin FS	
		%	42.1	33.7	50.5	4.20	8.40	Siemens/Dade Actin FSL	
		%	44.0	35.2	52.8	4.40	8.80	Siemens/Dade Pathromtin SL	
		%	44.0	35.2	52.8	4.40	8.80	Stago CK Prest	
Factor X activity	X因子活性	%	40.9	32.7	49.1	4.10	8.20	Stago Deficient VIII	
		%	40.3	32.2	48.4	4.05	8.10	Stago STA Factor VIII	
		%	75.8	60.6	91.0	7.60	15.2	HemosIL Factor X deficient plasma	
		%	72.4	57.9	86.9	7.25	14.5	HemosIL RecombiPlasTin 2G	
		%	86.4	69.1	104	8.80	17.6	Siemens Factor X Deficient Plasma	
		%	80.4	64.3	96.5	8.05	16.1	Siemens/Dade Innovin	
		%	88.8	71.0	107	9.10	18.2	Siemens/Dade Thromborel S	
		%	71.0	56.8	85.2	7.10	14.2	Stago Deficient Factor X	
		%	70.5	56.4	84.6	7.05	14.1	Stago STA-NeoPTimal	

Method 所有仪器均值		COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)					凝血质控品 (水平2)		
批号 (Lot. No): 405CG		货号 (Cat. No): CG5022		到期 (Expiry): 2025-09-28					
包装规格 (Size): 12 x 1 ml		Range			范围				
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	1SD	2SD	Method 方法	
Factor XI activity	XI因子活性	%	36.5	29.2	43.8	3.65	7.30	HemosIL Factor XI deficient plasma	
		%	35.5	28.4	42.6	3.55	7.10	HemosIL Synthasil	
		%	31.5	25.2	37.8	3.15	6.30	Siemens/Dade Actin FSL	
Factor XII activity	XII因子活性	%	68.5	54.8	82.2	6.85	13.7	HemosIL Factor XII deficient plasma	
		%	62.7	50.2	75.2	6.25	12.5	HemosIL Synthasil	
		%	68.4	54.7	82.1	6.85	13.7	Siemens Factor XII Deficient Plasma	
		%	68.7	55.0	82.4	6.85	13.7	Siemens/Dade Actin FS	
		%	69.4	55.5	83.3	6.95	13.9	Siemens/Dade Actin FSL	
		%	74.9	59.9	89.9	7.50	15.0	Siemens/Dade Pathromtin SL	
Fibrinogen	纤维蛋白原	g/l	2.04	1.63	2.45	0.205	0.410	Bio-Ksel Fibrinogen	
		mg/dl	204	163	245	20.5	41.0		
		g/l	2.01	1.61	2.41	0.200	0.400	Biolabo Bio-Fibri	
		mg/dl	201	161	241	20.0	40.0		
		g/l	1.96	1.57	2.35	0.195	0.390	Cobas Fibrinogen	
		mg/dl	196	157	235	19.5	39.0		
		g/l	1.61	1.29	1.93	0.160	0.320	Diagon Dia-F	
		mg/dl	161	129	193	16.0	32.0		
		g/l	1.69	1.35	2.03	0.170	0.340	DIRUI Fibrinogen	
		mg/dl	169	135	203	17.0	34.0		
		g/l	1.64	1.31	1.97	0.165	0.330	Erba Thrombin Reagent	
		mg/dl	164	131	197	16.5	33.0		
		g/l	2.06	1.65	2.47	0.205	0.410	HaemoDiagnostics Fibrinogen	
		mg/dl	206	165	247	20.5	41.0		
		g/l	2.13	1.70	2.56	0.215	0.430	Helena Clauss Fibrinogen 100	
		mg/dl	213	170	256	21.5	43.0		
		g/l	1.87	1.50	2.24	0.185	0.370	HemosIL Fibrinogen C	
		mg/dl	187	150	224	18.5	37.0		
		g/l	2.06	1.65	2.47	0.205	0.410	HemosIL PT-Fibrinogen	
		mg/dl	206	165	247	20.5	41.0		
		g/l	1.83	1.46	2.20	0.185	0.370	HemosIL PT-Fibrinogen HS Plus	
		mg/dl	183	146	220	18.5	37.0		
		g/l	1.77	1.42	2.12	0.175	0.350	HemosIL QFA (bovine thrombin)	
		mg/dl	177	142	212	17.5	35.0		
		g/l	1.33	1.06	1.60	0.135	0.270	HemosIL RecombiPlasTin 2G	
		mg/dl	133	106	160	13.5	27.0		
		g/l	2.06	1.65	2.47	0.205	0.410	Human Fibrinogen	
		mg/dl	206	165	247	20.5	41.0		
		g/l	1.98	1.58	2.38	0.200	0.400	Labitec Fibrinogen	
		mg/dl	198	158	238	20.0	40.0		

Method 所有仪器均值			COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)				凝血症质控品 (水平2)		
批号 (Lot. No): 405CG		货号 (Cat. No): CG5022	到期 (Expiry): 2025-09-28						
包装规格 (Size): 12 x 1 ml			Range		范围				
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	1SD	2SD	Method 方法	
Fibrinogen	纤维蛋白原	g/l	2.00	1.60	2.40	0.200	0.400	Media IVD Fibrinogen	
		mg/dl	200	160	240	20.0	40.0		
		g/l	1.90	1.52	2.28	0.190	0.380	Mindray/Longisland Fibrinogen	
		mg/dl	190	152	228	19.0	38.0		
		g/l	1.79	1.43	2.15	0.180	0.360	MTI Diagnostics Fibrinogen	
		mg/dl	179	143	215	18.0	36.0		
		g/l	1.62	1.30	1.94	0.160	0.320	Renam Fibrinogen	
		mg/dl	162	130	194	16.0	32.0		
		g/l	2.02	1.62	2.42	0.200	0.400	Roche Cobas t411 Fibrinogen	
		mg/dl	202	162	242	20.0	40.0		
		g/l	2.14	1.71	2.57	0.215	0.430	Roche Cobas t511/t711 Fibrinogen	
		mg/dl	214	171	257	21.5	43.0		
		g/l	2.11	1.69	2.53	0.210	0.420	Sclavo Fibrinogen	
		mg/dl	211	169	253	21.0	42.0		
		g/l	1.52	1.22	1.82	0.150	0.300	Siem/DadeThromboplastin deriv.Fibrinogen	
		mg/dl	152	122	182	15.0	30.0		
		g/l	1.48	1.18	1.78	0.150	0.300	Siemens/Dade Innovin derived Fibrinogen	
		mg/dl	148	118	178	15.0	30.0		
		g/l	1.91	1.53	2.29	0.190	0.380	Siemens/Dade Multifibren U	
		mg/dl	191	153	229	19.0	38.0		
		g/l	1.73	1.38	2.08	0.175	0.350	Siemens/Dade Thrombin	
		mg/dl	173	138	208	17.5	35.0		
		g/l	2.05	1.64	2.46	0.205	0.410	Stago Fibrinogen/ Liquid-Fib	
		mg/dl	205	164	246	20.5	41.0		
		g/l	1.99	1.59	2.39	0.200	0.400	Stago STA Fibri-Prest Automate	
		mg/dl	199	159	239	20.0	40.0		
		g/l	2.08	1.66	2.50	0.210	0.420	Tcoag TriniCLOT Fibrinogen	
		mg/dl	208	166	250	21.0	42.0		
		g/l	2.39	1.91	2.87	0.240	0.480	Technoclone Fibrinogen	
		mg/dl	239	191	287	24.0	48.0		
		g/l	1.80	1.44	2.16	0.180	0.360	Teco Teclot Fibrinogen	
		mg/dl	180	144	216	18.0	36.0		
		g/l	1.96	1.57	2.35	0.195	0.390	Thalys Fibrinogen	
		mg/dl	196	157	235	19.5	39.0		
		g/l	1.69	1.35	2.03	0.170	0.340	Yumizen G Fibrinogen	
		mg/dl	169	135	203	17.0	34.0		
Plasminogen activity	纤溶酶原活性	%	67.6	54.1	81.1	6.75	13.5	Hemosil Plasminogen	
		%	78.9	63.1	94.7	7.90	15.8	Siemens/Dade Berichrom Plasminogen	
Protein C activity	蛋白C活性	%	65.5	52.4	78.6	6.55	13.1	IL HemosIL Protein C	

Method 所有仪器均值		COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)					凝血质控品 (水平2)		
批号 (Lot. No): 405CG		货号 (Cat. No): CG5022		到期 (Expiry): 2025-09-28					
包装规格 (Size): 12 x 1 ml		Range			范围				
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	1SD	2SD	Method 方法	
Protein C activity	蛋白C活性	%	69.5	55.6	83.4	6.95	13.9	Siemens/Dade Berichrom Protein C	
		%	82.3	65.8	98.8	8.25	16.5	Siemens/Dade Protein C Clotting	
		%	74.6	59.7	89.5	7.45	14.9	STACHrom Protein C	
		%	74.8	59.8	89.8	7.50	15.0	Stago StaClot	
Protein S activity	蛋白S活性	%	71.2	49.8	92.6	10.7	21.4	IL HemosIL Free Protein S	
		%	76.5	53.6	99.4	11.5	22.9	IL HemosIL ProS	
		%	68.1	47.7	88.5	10.2	20.4	IL HemosIL Protein S Activity	
		%	70.9	49.6	92.2	10.7	21.3	Siemens Innovance Free PS Ag	
		%	57.1	40.0	74.2	8.55	17.1	Siemens/Dade Protein S Ac	
		%	75.1	52.6	97.6	11.3	22.5	STA Liatest	
		%	71.1	49.8	92.4	10.7	21.3	Stago StaClot	
Prothrombin Time (PT)	凝血酶原时间 (PT)	seconds	15.9	12.7	19.1	1.60	3.20	Beijing Succeeder PT	
		seconds	12.9	10.3	15.5	1.30	2.60	Bio-Ksel PT Plus	
		seconds	13.1	10.5	15.7	1.30	2.60	Bio-Ksel System PT	
		seconds	21.1	16.9	25.3	2.10	4.20	Biolabo Bio-TP	
		seconds	22.4	17.9	26.9	2.25	4.50	Biolabo Bio-TP (Low ISI)	
		seconds	19.7	15.8	23.6	1.95	3.90	Biomedica QuikCoag PT-HS	
		seconds	17.9	14.3	21.5	1.80	3.60	DG Lab APTT	
		seconds	17.9	14.3	21.5	1.80	3.60	Diagon, Dia-PT	
		seconds	15.8	12.6	19.0	1.60	3.20	Diagon, Dia-PT R	
		seconds	12.5	10.0	15.0	1.25	2.50	DIRUI PT	
		seconds	18.9	15.1	22.7	1.90	3.80	Dutch Diagnostics PT	
		seconds	18.1	14.5	21.7	1.80	3.60	Erba Protime PT	
		seconds	13.2	10.6	15.8	1.30	2.60	Fortress Diagnostics PT with Calcium	
		seconds	18.9	15.1	22.7	1.90	3.80	Helena Thromboplastin	
		seconds	20.2	16.2	24.2	2.00	4.00	HemosIL PT-Fibrinogen HS Plus	
		seconds	13.5	10.8	16.2	1.35	2.70	HemosIL ReadiPlasTin	
		seconds	13.6	10.9	16.3	1.35	2.70	HemosIL RecombiPlasTin 2G	
		seconds	18.6	14.9	22.3	1.85	3.70	Human Thromboplastin	
		seconds	17.9	14.3	21.5	1.80	3.60	Labitec PT	
		seconds	18.9	15.1	22.7	1.90	3.80	Media IVD PT	
		seconds	16.9	13.5	20.3	1.70	3.40	Mindray/Longisland PT	
		seconds	17.9	14.3	21.5	1.80	3.60	MTI Diagnostics Prothrombin Time	
		seconds	20.1	16.1	24.1	2.00	4.00	Roche Cobas PT Rec	
		seconds	14.3	11.4	17.2	1.45	2.90	Roche Cobas PT Screen	
		seconds	17.5	14.0	21.0	1.75	3.50	Sclavo PT	
		seconds	17.5	14.0	21.0	1.75	3.50	Siemens/Dade Innovin	
		seconds	13.8	11.0	16.6	1.40	2.80	Siemens/Dade Thromborel IS	
		seconds	13.6	10.9	16.3	1.35	2.70	Siemens/Dade Thromborel S	
		seconds	21.1	16.9	25.3	2.10	4.20	Stago Neoplastine CI	
		seconds	15.8	12.6	19.0	1.60	3.20	Stago Neoplastine R	
		seconds	19.8	15.8	23.8	2.00	4.00	Stago STA Neoplastine CI Plus	
		seconds	23.1	18.5	27.7	2.30	4.60	Stago STA-NeoPTimal	
		seconds	25.5	20.4	30.6	2.55	5.10	Tcoag TriniCLOT PT Excel S	

Method 所有仪器均值			COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)					凝血质控品 (水平2)	
批号 (Lot. No): 405CG		货号 (Cat. No): CG5022	到期 (Expiry): 2025-09-28						
包装规格 (Size): 12 x 1 ml			Range		范围				
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	1SD	2SD	Method 方法	
Prothrombin Time (PT)	凝血酶原时间 (PT)	seconds	19.8	15.8	23.8	2.00	4.00	Technoclone Technoplastin HIS	
		seconds	17.3	13.8	20.8	1.75	3.50	Teco TEClot PT-S	
		seconds	18.2	14.6	21.8	1.80	3.60	Thalys PT-HIS	
		seconds	17.4	13.9	20.9	1.75	3.50	Yumizen G PT	
		seconds	20.1	16.1	24.1	2.00	4.00	Yumizen G PT Liq	
Prothrombin Time (PT) activity	凝血酶原时间 (PT) 活性	%	85.8	68.6	103	8.60	17.2	Bio-Ksel PT Plus	
		%	49.3	39.4	59.2	4.95	9.90	Biolabo Bio-TP (Low ISI)	
		%	46.1	36.9	55.3	4.60	9.20	Diagon, Dia-PT	
		%	51.6	41.3	61.9	5.15	10.3	HemosIL PT-Fibrinogen HS Plus	
		%	75.5	60.4	90.6	7.55	15.1	HemosIL ReditiPlasTin	
		%	75.5	60.4	90.6	7.55	15.1	HemosIL RecombiPlasTin 2G	
		%	45.0	36.0	54.0	4.50	9.00	Human Thromboplastin	
		%	55.6	44.5	66.7	5.55	11.1	Media IVD PT	
		%	30.3	24.2	36.4	3.05	6.10	Roche Cobas PT Rec	
		%	65.4	52.3	78.5	6.55	13.1	Roche Cobas PT Screen	
		%	64.8	51.8	77.8	6.50	13.0	Sclavo PT	
		%	39.0	31.2	46.8	3.90	7.80	Siemens/Dade Innovin	
		%	72.7	58.2	87.2	7.25	14.5	Siemens/Dade Thromborel S	
		%	51.0	40.8	61.2	5.10	10.2	Stago Neoplastine CI	
		%	57.2	45.8	68.6	5.70	11.4	Stago Neoplastine CL Plus -Etalquick	
		%	72.4	57.9	86.9	7.25	14.5	Stago Neoplastine R	
		%	85.6	68.5	103	8.70	17.4	Stago SPA	
		%	52.7	42.2	63.2	5.25	10.5	Stago STA Neoplastine CI Plus	
		%	44.8	35.8	53.8	4.50	9.00	Stago STA-NeoPTimal	
		%	31.5	25.2	37.8	3.15	6.30	Tcoag TriniCLOT PT Excel S	
		%	71.6	57.3	85.9	7.15	14.3	Technoclone Technoplastin HIS	
		%	50.0	40.0	60.0	5.00	10.0	Teco TEClot PT-S	
		%	39.6	31.7	47.5	3.95	7.90	Thalys PT-HIS	
		%	46.6	37.3	55.9	4.65	9.30	Yumizen G PT	
		%	34.5	27.6	41.4	3.45	6.90	Yumizen G PT Liq	
Prothrombin Time (PT) INR	凝血酶原时间 (PT) 国际标准化比值	INR	1.27	1.02	1.52	0.125	0.250	Beijing Succeder PT	
		INR	1.19	0.952	1.43	0.120	0.240	Bio-Ksel PT Plus	
		INR	1.77	1.42	2.12	0.175	0.350	Biolabo Bio-TP	
		INR	1.83	1.46	2.20	0.185	0.370	Biolabo Bio-TP (Low ISI)	

Method 所有仪器均值			COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)					凝血质控品 (水平2)	
批号 (Lot. No): 405CG			货号 (Cat. No): CG5022		到期 (Expiry): 2025-09-28				
包装规格 (Size): 12 x 1 ml					Range	范围			
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	ISD	2SD	Method 方法	
Prothrombin Time (PT) INR	凝血酶原时间 (PT) 国际标准化比值	INR	1.70	1.36	2.04	0.170	0.340	Biomedica QuikCoag PT-HS	
		INR	1.54	1.23	1.85	0.155	0.310	DG Lab APTT	
		INR	1.49	1.19	1.79	0.150	0.300	Diagon, Dia-PT	
		INR	1.47	1.18	1.76	0.145	0.290	Diagon, Dia-PT R	
		INR	1.07	0.856	1.28	0.105	0.210	DIRUI PT	
		INR	1.62	1.30	1.94	0.160	0.320	Dutch Diagnostics PT	
		INR	1.52	1.22	1.82	0.150	0.300	Erba Prottime PT	
		INR	1.42	1.14	1.70	0.140	0.280	Helena Thromboplastin	
		INR	1.31	1.05	1.57	0.130	0.260	HemosIL PT-Fibrinogen	
		INR	1.51	1.21	1.81	0.150	0.300	HemosIL PT-Fibrinogen HS Plus	
		INR	1.21	0.968	1.45	0.120	0.240	HemosIL ReadiPlasTin	
		INR	1.20	0.960	1.44	0.120	0.240	HemosIL RecombiPlasTin 2G	
		INR	1.57	1.26	1.88	0.155	0.310	Human Thromboplastin	
		INR	1.47	1.18	1.76	0.145	0.290	Labitec PT	
		INR	1.45	1.16	1.74	0.145	0.290	Media IVD PT	
		INR	1.34	1.07	1.61	0.135	0.270	Mindray/Longisland PT	
		INR	1.66	1.33	1.99	0.165	0.330	MTI Diagnostics Prothrombin Time	
		INR	1.38	1.10	1.66	0.140	0.280	Pacific Thromboplastin DS	
		INR	2.08	1.66	2.50	0.210	0.420	Roche Cobas PT Rec	
		INR	1.37	1.10	1.64	0.135	0.270	Roche Cobas PT Screen	
		INR	1.42	1.14	1.70	0.140	0.280	Sclavo PT	
		INR	1.74	1.39	2.09	0.175	0.350	Siemens/Dade Innovin	
		INR	1.19	0.952	1.43	0.120	0.240	Siemens/Dade Thromborel S	
		INR	1.56	1.25	1.87	0.155	0.310	Spinreact PT	
		INR	1.80	1.44	2.16	0.180	0.360	Stago Neoplastine CI	
		INR	1.52	1.22	1.82	0.150	0.300	Stago Neoplastine CL Plus -Etaloquick	
		INR	1.17	0.936	1.40	0.115	0.230	Stago Neoplastine R	
		INR	1.07	0.856	1.28	0.105	0.210	Stago SPA	
		INR	1.61	1.29	1.93	0.160	0.320	Stago STA Neoplastine CI Plus	
		INR	1.80	1.44	2.16	0.180	0.360	Stago STA-NeoPTimal	
		INR	1.22	0.976	1.46	0.120	0.240	Technoclone Technoplastin HIS	
		INR	1.33	1.06	1.60	0.135	0.270	Teco TEClot PT-S	
		INR	1.60	1.28	1.92	0.160	0.320	Thalys PT-HIS	
		INR	3.20	2.56	3.84	0.320	0.640	Wiener Lab Soluplastin	
		INR	1.49	1.19	1.79	0.150	0.300	Yumizen G PT	
		INR	1.87	1.50	2.24	0.185	0.370	Yumizen G PT Liq	
		INR	1.73	1.38	2.08	0.175	0.350	Yumizen G PT Reco	
Prothrombin Time (PT) ratio	凝血酶原时间 (PT) 比值	RATIO	1.18	0.944	1.42	0.120	0.240	HemosIL RecombiPlasTin 2G	
		RATIO	1.53	1.22	1.84	0.155	0.310	Human Thromboplastin	

Method 所有仪器均值			COAGULATION CONTROL - LEVEL 2 (COAG CONTROL 2)					凝血质控品 (水平2)	
批号 (Lot. No): 405CG			货号 (Cat. No): CG5022		到期 (Expiry): 2025-09-28				
包装规格 (Size): 12 x 1 ml					Range	范围			
Analyte	分析项目	Unit 单位	Target 靶值	Low 低值	High 高值	1SD	2SD	Method 方法	
Prothrombin Time (PT) ratio	凝血酶原时间 (PT) 比值	RATIO	1.31	1.05	1.57	0.130	0.260	Sclavo PT	
		RATIO	1.77	1.42	2.12	0.175	0.350	Siemens/Dade Innovin	
		RATIO	1.18	0.944	1.42	0.120	0.240	Siemens/Dade Thromborel S	
		RATIO	1.58	1.26	1.90	0.160	0.320	Stago STA Neoplastine CI Plus	
		RATIO	1.76	1.41	2.11	0.175	0.350	Stago STA-NeoPTimal	
		RATIO	1.60	1.28	1.92	0.160	0.320	Thalys PT-HIS	
Thrombin Time (TT)	凝血酶时间 (TT)	seconds	25.9	20.7	31.1	2.60	5.20	Biolabo Bio-TT	
		seconds	19.6	15.7	23.5	1.95	3.90	Erba Thrombin Time	
		seconds	17.8	14.2	21.4	1.80	3.60	Helena Thrombin Clotting Time rgt.	
		seconds	23.6	18.9	28.3	2.35	4.70	HemosIL Thrombin Time 5ml	
		seconds	17.4	13.9	20.9	1.75	3.50	Mindray/Longisland TT	
		seconds	30.6	24.5	36.7	3.05	6.10	Renam TT	
		seconds	35.4	28.3	42.5	3.55	7.10	Siemens BC Thrombin	
		seconds	27.8	22.2	33.4	2.80	5.60	Siemens Test Thrombin	
		seconds	28.4	22.7	34.1	2.85	5.70	Siemens/Dade Thrombin	
		seconds	23.4	18.7	28.1	2.35	4.70	Siemens/Dade Thromboclotin	
		seconds	29.9	23.9	35.9	3.00	6.00	Stago Thrombin	
		seconds	27.1	21.7	32.5	2.70	5.40	Sysmex Thrombin Time	
		seconds	16.4	13.1	19.7	1.65	3.30	Thalys TT	