

适用细胞：

人骨髓间充质干细胞、人脂肪间充质干细胞、人脐带间充质干细胞、人羊膜间充质干细胞、人牙髓间充质干细胞、人脐带血间充质干细胞，其他组织来源的人间充质干细胞

产品特点：

操作简便，稳定性高，服务型产品

产品：人间充质干细胞成骨诱导试剂盒

货号：FY200006

诱导周期：1-3周

自带判定方法：茜素红染色

产品：人间充质干细胞成脂诱导试剂盒

货号：FY200007

诱导周期：1-3周

自带判定方法：油红O染色

产品：人间充质干细胞成软骨诱导试剂盒

货号：FY200008

诱导周期：2-3周

自带判定方法：阿利新蓝染色



扫二维码
关注更多

弗元（上海）生物科技有限公司




网址：www.fuyuanbio.com

邮箱：info@fuyuanbio.com

QQ：3376176889、842394155、262017109

电话：021-80370389

引用实例：



引用文献：

Yue Hong, Shan Shan, Ye Gu, Haidong Huang, Quncheng Zhang, Yang Han, Yongpin Dong, Zeyu Liu, Moli Huang, Tao Ren. Malfunction of airway basal stem cells plays a crucial role in pathophysiology of tracheobronchopathia osteoplastica. *Nat Commun.* 2022 Mar 14;13(1):1309. (IF 14.919)

Li, T., Xu, Y., Wang, Y. et al. Differential expression profiles of long noncoding RNAs and mRNAs in human bone marrow mesenchymal stem cells after exposure to a high dosage of dexamethasone. *Stem Cell Res Ther* 12, 9 (2021). (IF 8.079)

Zhang L, Zhang C, Ji Y, et al. Effects of Z-value on physicochemical and biological properties of β -SiAlONs ceramics[J]. *Ceramics International*, 2021. (IF 5.532)

Sun J, Zhang F, Luo X, Shi G, Li F, Zheng B, Guo Y, Shi J, Li L. Long noncoding RNA AC092155 facilitates osteogenic differentiation of adipose-derived stem cells through the miR-143-3p/STMN1 axis. *J Gene Med.* 2021 Aug;23(8):e3363. (IF 4.152)

Zhao ZH, Ma XL, Zhao B, Tian P, Ma JX, Kang JY, Zhang Y, Guo Y, Sun L. Naringin-inlaid silk fibroin/hydroxyapatite scaffold enhances human umbilical cord-derived mesenchymal stem cell-based bone regeneration. *Cell Prolif.* 2021 Jul;54(7):e13043. (IF 8.755)



扫二维码
关注更多

弗元（上海）生物科技有限公司

网址：www.fuyuanbio.com

邮箱：info@fuyuanbio.com

QQ：3376176889、842394155、262017109

电话：021-80370389